Paediatric positioning system



The multifunctional postural solution for kids



Destination of use

Ideal for kids and children affected by inadeguate postural control as results of several pathologies. Ermes allows customisation and adjustments during the assessment, delivery and after-delivery follow-up, for all users affected by leg length discrepancies, wind swept, pelvic obliquity, anterior/posterior pelvis tilt even combined and/or related to scoliosis, lordosis and kyphosis.



Integrated Positioning System ERMES + TIPCO



Main Features ERMES + TIPCO



Integrated Positioning System ERMES + ADACTA GIÒ



Main Features ERMES + ADACTA GIÒ



Ermes Materials

• Seat and Back positioning system with aluminium shell and polyethylene inserts kit.

• Two layers of foam for the backrest (self-modeling polymeric foam and reticulated foam for ventilation) and microclimatic ventilated cover.

- Kit of covers for the backrest in microclimatic material, elastic and breathable.
- Seat padding in visco-elastic foam.
- FREE SHAPE microclimatic and incontinent padded cover kit for the seat.

Guide to the choice the right Seat



Optional Parts and Accessories (with or without up-charge. Please refer to the order form)

For partially sighted users this document can be visualized in PDF format on the web site www.promedicare.eu with the possibility to enlarge the text

Technical Information

	ERMES		Size 2	Size 3	Size 4
	SEAT WIDTH (cm)		≤25	18-28	22-35
	SEAT DEPTH (cm)		19-26	22-34	31-43
	TRUNK HEIGHT (cm)		26-38	33-45	41-53
	CHEST WIDTH (cm)	BRACKET A	19-25	22-28	27-33
		BRACKET B	14-20	17-23	22-28
	LEG MAXIMUM ABDUCTION (°)		16	16	16
	LEG MAXIMUM ADDUCTION (°)		16	16	16
	HEIGHT OF POPLITEAL SPACE TO FOOTREST (cm)		8-19	8-29	8-39
	ARM SUPPORT HEIGHT (cm)		14-19	20-26	20-26
	BACKREST RECLINING (°)		-3/+20	-3/+20	-3/+20
	ERMES LOAD (kg)		25	35	50
	ERMES WEIGHT (kg)		5.2	8.1	12.1
Oę	ERMES TIPCO		Size 2	Size 3	Size 4
	MAXIMUM SEAT TO FLOOR HEIGHT (cm)		63	63	63
	WIDTH OVERALL DIMENSIONS (cm)		67	67	67
	DEPTH OVERALL DIMEN	DEPTH OVERALL DIMENSIONS (cm)		80	88
ETTERA DI POSTURA II	LOAD (kg)			35	50
	TIPCO WEIGHT (kg) TOTAL WEIGHT (kg)		13.6	13.6	14
$\mathbf{\Sigma}$			18.8	21.7	26.1
	TILT IN SPACE (°)		-5/+30	-5/+30	-5/+30
					1
	ERMES GIÒ		Size 2	Size 3	Size 4
	SEAT TO FLOOR HEIG	GHT (cm)	52	52	53
	MAX WIDTH OVERALL DIMENSIONS (WI	TH 400/500 mm WHEELS) (cm)	48	53	58
Ŭ Ŭ	MIN/ MAX DEPTH OVERALL D	IMENSIONS (cm)	73-100	75-102	81-109
	LOAD (kg)		25	35	50

MAX WIDTH OVERALL DIMENSIONS (WITH 400/500 mm WHEELS) (cm)		53	58
MIN/ MAX DEPTH OVERALL DIMENSIONS (cm)		75-102	81-109
LOAD (kg)		35	50
GIÒ WEIGHT (STANDARD CONFIGURATION) (kg)		13.4	14.5
TOTAL WEIGHT (kg)	18	21.5	26.6
PREDETERMINED BACK POSTS RECLINING (°)	3 - 10 - 18		
TILT IN SPACE (°) *Upgraded tilt in space of 4° with an increment of the starting position of the frame by operating on the positioning holes of the caster fork	2-35 (da 6 a 39)*		

Manufactured by:

edicare

Via Montagna Z.I. Lotto 41 72023 Mesagne (Br) - ITALY Tel. +39 0831 777840 - Fax. +39 0831 730739 Website: www.promedicare.eu e-mail: info@promedicare.it

> CERTIFIED COMPANY WITH QUALITY MANAGEMENT SYSTEM ISO 9001/ISO 13485





The contents of this catalogue is the property of Pro Medicare s.r.l. Any improper or anauthorised use is strictly forbidden www.promedicare.eu

The manufacturer reserves the right to make changes without notice.