

ORDER FORM

Invoicing Details Company: Address:.....N°:.....P.Office: Town: State: Tel: Fax: V.A.T.#: Order #: Date of the order...../...../.....	Shipping Details Ship to : <p style="text-align: center;">CUSTOMER STAMP AND SIGNATURE</p>
--	---

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Technical Standard Features

- Frame Adacta Liga Light anodised satin
- Folding down Back posts
- Reclining mechanism by predetermined plate w/pin 3°,9°,15°,22°
- Continuous Tilt in Space by drive lever from 2° up to 35° (from 6° to 39°)¹
- ¹ Tilt in space can be up-graded of 4° with an increment of the starting position of the frame by operating on the positioning holes of the caster fork
- Adjustment of the seat depth
- Height adjustment: High (49 cm from the floor) / low (43 cm from the floor)*
- * (with 300mm and 500mm wheels only)
- Footplates with 90° fixed knee angle and detachable hangers
- Anti-tip system fixed
- Parking brake with lever
- Rigid wooden seat base covered leather flex
- Multi-Adjustable padded hip guides
- Multi-Adjustable padded arm supports
- Versa Inserto Trunk/Pelvis positioning system
- Capitis Headrest
- Fixatis Pelvic Belt



The images are purely indicative and may not reflect the standard configuration.



◊ = Standard ○ = Optional N/C □ = Optional with up-charge

* If a selection is not made, the product will be delivered in its standard configuration.

Frame Colours	Frame Type	Frame Sizes (Outside tubes Width x Depth cm.)
◊ Turquoise ○ Pink	◊ Rigid ○ Folding	○ US(30x38) ○ XXS(35x38) ○ XS(35x42) ○ XS0(38x42)
Knee Angle	◊ 90° ○ 75° □ elevating	
Castors	◊ 150 mm (6") poly grey ○ black ○ 175 mm (7") poly ○ grey ○ black	<i>(not suitable with 600mm (24") posterior wheels)</i>
Posterior Wheels (height 43 cm from seat to the floor with 300 mm and 500 mm wheels; height 47 cm from seat to the floor with 600 mm wheels) <i>Note: If not indicated, the wheels will be gray</i>		
◊ 300 mm (12") poly grey ○ black ○ 300 mm (12") pneumatic		
□ 300 mm (12") pneumatic w/drum brake + lever □ grey □ black		
○ 500 mm (20") pneumatic without handrim □ handrims for 500 mm wheels		
□ 500 mm (20") pneumatic without handrim and w/drum brake + lever		
□ 600 mm (24") pneumatic w/handrim (not available on size US) (Note: The tilt in space is reduced to 20°, height from the floor 47 cm)		
Anti-tip system detachable		□
Pelvis Positioning System		
◊ Inserto Seat Basico	□ Inserto Seat Mediano	□ Inserto Seat Domino
Trunk Positioning System		
◊ Inserto Back Effecto Lato	○ Inserto Back Supporto	
Backposts (Note: Push handles height 63 cm.)		

- push handles height adjustable by pin (max adjustment 10 cm with 2,5 cm gap)
- telescopic push handles kit with lever block (max adjustment 10 cm)
- angle-adjustable reclining single handlebar kit for rigid frame
- angle-adjustable reclining handlebar kit for rigid frame

Accessories

Thoracic Supports for Inserto Back Supporto *cross the type, the shapes and the size desired

Type of Bracket	Shape of the Pad		Dimensions (cm)				
Fixed <input type="checkbox"/> RH <input type="checkbox"/> LH	Flat <input type="checkbox"/> RH <input type="checkbox"/> LH		cm	<input type="checkbox"/> size 1	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4
Swing-Away <input type="checkbox"/> RH <input type="checkbox"/> LH	Contoured <input type="checkbox"/> RH <input type="checkbox"/> LH		H	6	8	10	12
			L	6,5	8,5	10,5	13,5
			cm	<input type="checkbox"/> size 1	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	
			H	6	8	10	
			L	10,5	12,5	14,5	

Thoracic Supports for Inserto Back Effecto Lato *cross the type, the shape of the brackets, the pad and the size desired

Type of Bracket	Shape of the Pad		Dimensions (cm)				
Fixed <input type="checkbox"/> RH <input type="checkbox"/> LH	Flat <input type="checkbox"/> RH <input type="checkbox"/> LH		cm	<input type="checkbox"/> size 1	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	<input type="checkbox"/> size 4
Swing-Away <input type="checkbox"/> RH <input type="checkbox"/> LH	Contoured <input type="checkbox"/> RH <input type="checkbox"/> LH		H	6	8	10	12
Detachable <input type="checkbox"/> RH <input type="checkbox"/> LH	Shape of the brackets *for a right selection see TECHNICAL INFORMATION Type A: <input type="checkbox"/> RH <input type="checkbox"/> LH Type B: <input type="checkbox"/> RH <input type="checkbox"/> LH		L	8	10	12	15
			cm	<input type="checkbox"/> size 1	<input type="checkbox"/> size 2	<input type="checkbox"/> size 3	
			H	6	8	10	
			L	12	14	16	

Capitis Headrest Comfort <input type="radio"/> Small <input checked="" type="radio"/> Medium <input type="radio"/> Large	<input type="checkbox"/> Kit of ABS laterals height adjustable w/armrest
Fixatis Pelvic belt 2P.2A. with squeeze plastic buckle size: <input type="radio"/> XS <input checked="" type="radio"/> S <input type="radio"/> M	Fixatis Dynamic Butterfly Harness <input type="checkbox"/> XS <input type="checkbox"/> S <input type="checkbox"/> M
<input type="checkbox"/> Separate Calf Supports swing-away	<input type="checkbox"/> Padded Calf Strap
<input type="checkbox"/> Detachable Abduction Block	<input type="checkbox"/> Pedal for barriers <input type="checkbox"/> RH <input type="checkbox"/> LH
<input type="checkbox"/> Rectangular Wooden Tray w/lateral flip hardware <input type="checkbox"/> Rectangular Wooden Tray w/detachable hardware <input type="checkbox"/> Padding for Tray	<input type="checkbox"/> Vehicle Transportation Kit <i>Only with 300mm and 500mm posterior wheels Do not use the wheelchair as a seat in the vehicle if the weight of the user is lower than 22 kg.</i>

Sizes of the arm supports and hip guides. **NOTE: The components cannot be interchangeable**

Linear Hip Guides						ABS Lateral					
	cm	US	XXS	XS	XSO		cm	US	XXS	XS	XSO
	H	10	10	10	10		H	11	11	11	13
	L	25	25	30	30		L	24	24	24	28
Flat Arm Support						Armrest of the ABS Lateral kit					
	cm	US	XXS	XS	XSO		cm	US	XXS	XS	XSO
	H	10	10	10	10		H	5	5	5	5
	L	25	25	30	30		L	27	27	27	33

Technical Information

Titti Size	Maximum Load (kg)	Weight of the Frame (kg)	Pivot Width (cm)	Maximum overall dimensions in width / length (90° knee angle and footrest positioned forwards) Values expressed in cm.					
				Wheel 300 strd.	Wheel 300 w/drum brake.	Wheel 500 strd.	Wheel 500 w/drum brake.	Wheel 600	
US	Rigid	50	14.1	69	48/96	49/96	50/101	53/101	/
	Folding	50	14.9						
XXS	Rigid	50	14.4	69	53/96	54/96	55/101	58/101	55/104
	Folding	50	15.3						
XS	Rigid	50	14.7	72	53/99	54/99	55/102	58/102	55/107
	Folding	50	15.6						
XSO	Rigid	50	15.4	75	56/99	57/99	58/102	61/102	58/107
	Folding	50	15.9						

Note: Frame weight only with wheels, knee angle and footrest as standard configuration. The length dimensions, referred to the average values, are subjected to variations due to the displacement of the wheel plates. The values reported are intended for the standard configuration.

Titti Size	User Parameters		Frame Depth Adjustments (cm)	Cushion Size (W x D) (cm)	Supporto Back sizes (cm)				Effecto Lato Back sizes (cm)			
	Pelvis Width (cm)	Seat Depth (cm)			Backrest	Height	Shell Width	Trunk Width	Backrest	Height	Shell Width	Trunk Width
US	20-30	28-38	36 to 44	30x38	3033	33	27	17-23	3033	33	30	w/bracket A 18-28 w/bracket B 13-23
XXS	26-36	28-38	36 to 44	36x38	3336	36	30	20-26	3338	38	33	w/bracket A 21-31 w/bracket B 16-26
XS	26-36	32-42	40 to 48	36x42	3640	40	33	23-29	3642	42	36	w/bracket A 24-34 w/bracket B 19-29
XSO	29-39	34-44	40 to 48	39x45	3944	44	36	26-34	3945	45	39	w/bracket A 27-37 w/bracket B 22-32